

## Choose Your Bedpan Washer

### Top Loading Bedpan Washers

Features	Hands Free	Antimicrobial Lid	Capacity	Cycle Time *	Water Supply	Rim Flush	Scale Mgmt System
Mini	✓	✓	1 Bedpan 1 Urine Bottle	6 mins	Cold	✗	✓
Midi	✓	✓	2 Bedpans 1 Urine Bottle	7 mins	Cold	✗	✓
Optima 2	✓	✓	2 Bedpans 2 Urine Bottles	8 mins	Hot & Cold	✗	✓
Maxi	✓	✗	3 Bedpans 2 Urine Bottles	7 mins	Cold	✗	✓
Maxi +	✓	✗	3 Bedpans 2 Urine Bottles	7 mins	Cold	✓	✓
Optima 3	✓	✗	3 Bedpans 3 Urine Bottles	10 mins	Hot & Cold	✗	✓

### Front loading Bedpan Washers

Features	Hands Free	Antimicrobial Lid	Capacity	Cycle Time *	Water Supply	Rim Flush	Scale Mgmt System
XLC	✓	✗	1 Bedpan 3 Urine Bottles	7-8 mins	Hot & Cold	✗	✓
XLS	✓	✗	1 Bedpan 3 Urine Bottles	7-8 mins	Hot & Cold	✗	✓
XL2	✓	✗	2 Bedpans 2 Urine Bottles	7-8 mins	Hot & Cold	✗	✓

\*Cycle times may vary depending on water pressure

### Why Choose a DDC Dolphin Bedpan Washer Disinfector?



Automatic scale management system, suitable for soft and hard water areas, to prevent blockages.



Optional programmable audio in your own language to ensure correct and proper use by staff.



Antimicrobial technology embedded in lid, stopping the growth of harmful bacteria in its tracks.



Completely hands-free, no touch operation, reducing the opportunity for infections to spread to clinicians and further patients.

Get in touch with one of our washer disinfector experts.

+44 (0) 01202 731555 | info@ddcdolphin.com | www.ddcdolphin.com



## Bedpan Washer Disinfectors

### Redefining Infection Control

Our range of washer disinfectors have been developed over decades to provide you with a fast and safe human waste disposal system with passive safety measures to minimise the risk of infection. Designed to empty, wash and disinfect human waste containers such as bedpans, commode pots and urine bottles, our bedpan washer disinfectors incorporate the latest technological advances and comply with HTM2030 and HTM0101 regulations. Our range of both top and front loading machines provide solutions for all your requirements.

### Benefits of DDC Dolphin Bedpan Washer Disinfectors

- 1 Reduce cross contamination risk with hands free technology.
- 2 Large, medium and smaller capacity models to suit various environments.
- 3 Automatic thermal disinfection to HTM2030 standards.
- 4 Extremely thorough wash performance.



# Redefining Infection Control

With superior and considerate build quality DDC Dolphin are raising the standards in the fight against infection.

We use only the highest grade components and materials meaning you are safe in the knowledge that your machine has been built to last.



## Patented antimicrobial surfaces

– to minimise cross infection and reduce the possibility of microbial growth.



## High grade stainless steel casing

– which is robust and long-lasting, not easily corroded, and easy to clean and maintain.



## Standard electrical supply

– so you can plug in and go!



## Short washing disinfection cycle

– eco-friendly and cost effective.



**Compact** – to save space, making it suitable for all room types and sizes.



## Patented hands-free open & close technology

– using a foot cup or infrared sensor to minimise cross infection and make it easy to use so clinicians can quickly return to caring for patients.

**Top or front loading** – for a simple loading and emptying process.



## Range and capacity

– can cater for any workload or load configuration. Our Maxi machine has the largest capacity in the industry.



**Unique multi-utensil racking design** – extremely convenient with one rack to fit all products.



**Scale Management** – automatic scale management system, suitable for soft and hard water areas.

Get in touch with one of our washer disinfector experts.



+44 (0) 01202 731555



info@ddcdolphin.com



www.ddcdolphin.com



**DDC**  
DOLPHIN